

# ENDOCUFF® Compatibility Schedule 2015

FUJINON		
Scope Model		Vision Model
EC-600WM/WI/WL	New series	AEC110 (Blue)
EC-590WM4/WI4/WL4/WI/WLY	Adult	AEC110 (Blue)
EC-590ZW3 M/L	Adult Zoom	AEC120 (Green)
EC-530FI/FL	Adult	AEC120 (Green)
EC-530WM3/WI3/WL3/WI	Adult	AEC120 (Green)
EC-450 Adult Series	Adult	AEC120 (Green)
EC-200MR	Adult	AEC120 (Green)
EC-250 Adult Series	Adult	AEC120 (Green)
EC-201WI	Adult	AEC120 (Green)
EC-530 DM/DL	Adult	AEC120 (Green)
EC-600 ZW/M/L	Adult	AEC120 (Green)
EC-530MP/LP	Paediatric	AEC130 (Purple)
EC-530LS	Intermediate	AEC130 (Purple)
EC-450 Paediatric Series	Paediatric	AEC130 (Purple)
EC-450LP/LS5	Intermediate	AEC130 (Purple)
EC-250 Intermediate Series	Intermediate	AEC130 (Purple)
OLYMPUS		
Scope Model		Vision Model
CF-H290L/I	Adult HD	AEC110 (Blue)
CF-Q260AI/AL/DI/DL	Adult Scope Guide	AEC110 (Blue)
CF-240AI/AL/DI/DL/I/L	Adult	AEC110 (Blue)
PCF-H190DL/I	Paediatric	AEC110 (Blue)
CF-H185L/I	Adult	AEC110 (Blue)
CF-H/Q170L/I	Adult	AEC110 (Blue)
CF-145 Series	Adult	AEC110 (Blue)
CF-HQ290L/I	Adult HD Dual Focus	AEC120 (Green)
CF-H260AI/AL/AZI/AZL	Adult HD Zoom	AEC120 (Green)
CF-Q240AI/AL/I/L/ZI/ZL	Adult	AEC120 (Green)
CF-HQ190L/I	Adult HD	AEC120 (Green)
CF-H190L/I	Adult HD	AEC120 (Green)
CF-Q180 AL/I	Adult	AEC120 (Green)
CF-H180 AL/I	Adult HD	AEC120 (Green)
CF-H260DL/I	Adult HD Scope Guide	AEC120 (Green)
CF-165 Series	Adult	AEC120 (Green)
CF-160 Series	Adult	AEC120 (Green)
CF-Q140 Series	Adult	AEC120 (Green)
CF-1T140L	Adult	AEC120 (Green)
PCF-H290L/I	Paediatric HD	AEC130 (Purple)
PCF-Q260AZI	Paediatric Zoom	AEC130 (Purple)
PCF-Q260AL/I	Paediatric	AEC130 (Purple)
PCF-H190L/I	Paediatric	AEC130 (Purple)
PCF-H180AL/I	Paediatric HD	AEC130 (Purple)
PCF-Q180	Paediatric	AEC130 (Purple)
PCF-160 Series	Paediatric	AEC130 (Purple)
CF-H180DL/I	Adult HD (Scope Guide)	AEC140 (Orange)
PCF-140L	Paediatric	Not Compatible
PCF/PH190L/I	Ultraslim	Not compatible
PCF-PQ260L/I	Adult	Not compatible
PENTAX		
Scope Model		Vision Model
EC38-i10M/M2/F/F2/L	Adult	AEC110 (Blue)
EC-3890TMK/MK2/FK/FK2/LK	Adult	AEC110 (Blue)

EC-3890MZI/FZI/LZI/SI2/TFK	Adult	AEC110 (Blue)
EC-3970LK	Adult	AEC120 (Green)
EC-3890FI/LI/Fi2(HD)	Adult	AEC120 (Green)
EC-3890LK	Adult	AEC120 (Green)
EC-3890TLK	Dual channel Adult	AEC120 (Green)
Scope Model		Vision Model
EC-3880FK2	Adult	AEC120 (Green)
EC-3870LK/LZK/TLK/FK2/MK2	Adult	AEC120 (Green)
EC-3840FK2	Adult	AEC120 (Green)
EC-3830LK/P/FK2	Adult	AEC120 (Green)
EC-380LKP/FKP/FK2P/MK2P	Adult	AEC120 (Green)
EC-38LV	Adult	AEC120 (Green)
EC-3940FK	Paediatric	AEC130 (Purple)
EC-3490 FI	Paediatric	AEC130 (Purple)
EC-3490 MK/FK/LK/LI/TLI/TFI	Paediatric	AEC130 (Purple)
EC-3470LK	Paediatric	AEC130 (Purple)
EC-34-i10M/F/L	Paediatric	AEC130 (Purple)
EC-3490TLF/MI	Retroview	Not compatible
EC-2990LI	Ultraslim	Not compatible
Scope Model		Vision Model
<b>STORZ</b>		
13924-PKS	Adult	AEC120 (Green)
13925-PKS	Adult	AEC120 (Green)
13907-PKS	Paediatric	AEC130 (Purple)
13906-PKS	Paediatric	AEC130 (Purple)
13905-PKS	Adult	AEC140 (Orange)
13904-PKS	Adult	AEC140 (Orange)

RELEASE: NOV 2015