


# ENDOCUFF VISION<sup>®</sup>



The hiding is over...



## ENDOCUFF VISION®

**ENDOCUFF VISION® has been designed to give an enhanced view of the entire colon. A single row of longer arms gently evert and flatten folds resulting in statistically significant and clinically relevant improvements in adenoma detection rates, as compared with unassisted colonoscopies.**

ENDOCUFF VISION® is disposable and simply pushes over the tip of all commonly used colonoscopes. Its unique hinged arms fall flat against the shaft of the colonoscope to create a smooth low friction surface during forward advancement. During withdrawal, the arms flare out and stabilize the tip, gently stretching the mucosal surface for close inspection.

### ENDOCUFF VISION® :

- Delivers yet more tip control throughout the colon without compromising intubation and improving loop management.
- Gives an early and controlled view of the upstream surface of large folds so no need for repeated intubation.
- Enables detailed inspection of flexures by preventing sudden slip back and red out.
- Opens concertinaed tracts of sigmoid colon as the arms gently separate the folds to provide a clear view of the mucosa between them.
- Acts as a handbrake, optimising tip position during therapy and polyp retrieval.
- There really is no hiding place for polyps.