

## ENDOCUFF VISION® Compatibility Schedule June 2017

Scope Model	FUJINON		Vision Model
EC-760R-V/L	Adult		ARV110 (Blue)
EC-760ZP-VM/VL	Adult		ARV110 (Blue)
EC-600WM/WI/WL	New series		ARV110 (Blue)
EC-590WM4/WI4/WL4/WI/WLY	Adult		ARV110 (Blue)
EC-590ZW3 M/L	Adult Zoom		ARV120 (Green)
EC-530FI/FL/FM	Adult		ARV120 (Green)
EC-530WM3/WI3/WL3/WI	Adult		ARV120 (Green)
EC-450 Adult Series	Adult		ARV120 (Green)
EC-200MR	Adult		ARV120 (Green)
EC-250 Adult Series	Adult		ARV120 (Green)
EC-201WI	Adult		ARV120 (Green)
EC-530 MT/IT/LT	Adult		ARV120 (Green)
EC-530 DM/DL	Adult		ARV120 (Green)
EC-600 ZW/M/L	Adult		ARV120 (Green)
EC-530MP/LP	Paediatric		ARV130 (Purple)
EC-530LS	Intermediate		ARV130 (Purple)
EC-450 Paediatric Series	Paediatric		ARV130 (Purple)
EC-450LP/LS5	Intermediate		ARV130 (Purple)
EC-250 Intermediate Series	Intermediate		ARV130 (Purple)
Scope Model	OLYMPUS		Vision Model
CF-H290L/I	Adult HD		ARV110 (Blue)
CF-Q260AI/AL/DI/DL	Adult Scope Guide		ARV110 (Blue)
CF-240AI/AL/DI/DL/I/L	Adult		ARV110 (Blue)
PCF-H190DL/I	Paediatric		ARV110 (Blue)
CF-H185L/I	Adult		ARV110 (Blue)
CF-H/Q170L/I	Adult		ARV110 (Blue)
CF-145 Series	Adult		ARV110 (Blue)
CF-HQ290L/I	Adult HD Dual Focus		ARV120 (Green)
CF-H260AI/AL/AZI/AZL	Adult HD Zoom		ARV120 (Green)
CF-Q240AI/AL/I/L/ZI/ZL	Adult		ARV120 (Green)
CF-HQ190L/I	Adult HD		ARV120 (Green)
CF-H190L/I	Adult HD		ARV120 (Green)
CF-Q180 AL/I	Adult		ARV120 (Green)
CF-H180 AL/I	Adult HD		ARV120 (Green)
CF-H260DL/I	Adult HD Scope Guide		ARV120 (Green)
CF-165 Series	Adult		ARV120 (Green)
CF-160 Series	Adult		ARV120 (Green)
CF-Q140 Series	Adult		ARV120 (Green)
CF-1T140L	Adult		ARV120 (Green)

**Note: For Canada Only**

Not all medical devices listed are licenced in accordance with Canadian law. Please contact your Olympus representative to find out what is registered for sale in your market.

OLYMPUS		
Scope Model		Vision Model
PCF-H290L/I	Paediatric HD	ARV130 (Purple)
PCF-Q260AZI	Paediatric Zoom	ARV130 (Purple)
PCF-Q260AL/I	Paediatric	ARV130 (Purple)
PCF-H190L/I	Paediatric	ARV130 (Purple)
PCF-H180AL/I	Paediatric HD	ARV130 (Purple)
PCF-Q180	Paediatric	ARV130 (Purple)
PCF-160 Series	Paediatric	ARV130 (Purple)
CF-H180DL/I	Adult HD (Scope Guide)	ARV140 (Orange)
PCF-140L	Paediatric	Not Compatible
PCF/PH190L/I	Ultraslim	Not compatible
PCF-PQ260L/I	Adult	Not compatible
PENTAX		
Scope Model		Vision Model
EC38-i10M/M2/F/F2/L/NF	Adult	ARV110 (Blue)
EC-3890TMK/MK2/FK/FK2/LK	Adult	ARV110 (Blue)
EC-3890MZI/FZI/LZI/SI2/TFK	Adult	ARV110 (Blue)
EC-3890FI/LI/Fi2(HD)/Mi2	Adult	ARV110 (Blue)
EC-3890TLK	Dual Channel Adult	ARV110 (Blue)
EC-3970LK	Adult	ARV120 (Green)
EC-3880FK2	Adult	ARV120 (Green)
EC-3870LK/LZK/TLK/FK2/MK2	Adult	ARV120 (Green)
EC-3840FK2	Adult	ARV120 (Green)
EC-3830LK/P/FK2	Adult	ARV120 (Green)
EC-380LKP/FKP/FK2P/MK2P	Adult	ARV120 (Green)
EC-38LV	Adult	ARV120 (Green)
EC-3940FK	Paediatric	ARV130 (Purple)
EC-3490 FI	Paediatric	ARV130 (Purple)
EC-3490 MK/FK/LK/LI/TLI/TFI	Paediatric	ARV130 (Purple)
EC-3470LK	Paediatric	ARV130 (Purple)
EC-34-i10M/F/L	Paediatric	ARV130 (Purple)
EC-3490TLF/MI	Retroview	Not compatible
EC-2990LI	Ultraslim	Not compatible
STORZ		
Scope Model		Vision Model
13924-PKS	Adult	ARV120 (Green)
13925-PKS	Adult	ARV120 (Green)
13907-PKS	Paediatric	ARV130 (Purple)
13906-PKS	Paediatric	ARV130 (Purple)
13905-PKS	Adult	ARV140 (Orange)
13904-PKS	Adult	ARV140 (Orange)

**Note: For Canada Only**

Not all medical devices listed are licenced in accordance with Canadian law. Please contact your Olympus representative to find out what is registered for sale in your market.