

SCOPE MODEL	FUJINON	ENDOCUFF VISION MODEL
EC-590WM4/WI4/WL4/WI/WLY	ADULT	ARV110 (BLUE)
EC-600WM/WI/WL	NEW SERIES	ARV110 (BLUE)
EC-760R-V/L	ADULT	ARV110 (BLUE)
EC-760ZP-VM/VL	ADULT	ARV110 (BLUE)
EC-200MR	ADULT	ARV120 (GREEN)
EC-201WI	ADULT	ARV120 (GREEN)
EC-250	ADULT SERIES	ARV120 (GREEN)
EC-450	ADULT SERIES	ARV120 (GREEN)
EC-530 DM/DL	ADULT	ARV120 (GREEN)
EC-530 MT/IT/LT	ADULT	ARV120 (GREEN)
EC-530 FI/FL/FM	ADULT	ARV120 (GREEN)
EC-530WM3/WI3/WL3/WI	ADULT	ARV120 (GREEN)
EC-590ZW3 M/L	ADULT ZOOM	ARV120 (GREEN)
EC-600 ZW/M/L	ADULT	ARV120 (GREEN)
EC-250	INTERMEDIATE SERIES	ARV130 (PURPLE)
EC-450	PAEDIATRIC SERIES	ARV130 (PURPLE)
EC-450LP/LS5	INTERMEDIATE	ARV130 (PURPLE)
EC-530LS	INTERMEDIATE	ARV130 (PURPLE)
EC-530MP/LP	PAEDIATRIC	ARV130 (PURPLE)
SCOPE MODEL	OLYMPUS	ENDOCUFF VISION MODEL
CF-240AL/I	ADULT	ARV110 (BLUE)
CF-240DL/I	ADULT	ARV110 (BLUE)
CF-240L/I	ADULT	ARV110 (BLUE)
CF-H185L/I	ADULT	ARV110 (BLUE)
CF-H290L/I	ADULT HD	ARV110 (BLUE)
CF-Q145L/I	ADULT	ARV110 (BLUE)
CF-Q260AL/I	ADULT SCOPE GUIDE	ARV110 (BLUE)
CF-Q260DL/I	ADULT SCOPE GUIDE	ARV110 (BLUE)
PCF-H190DL/I	PAEDIATRIC	ARV110 (BLUE)
CF-140L/I	ADULT	ARV120 (GREEN)
CF-140S	ADULT	ARV120 (GREEN)
CF-1T140L/I	ADULT	ARV120 (GREEN)
CF-H170L/1	ADULT	ARV120 (GREEN)
CF-H180 AL/I	ADULT HD	ARV120 (GREEN)
CF-H190L/I	ADULT HD	ARV120 (GREEN)
CF-H260AL/I	ADULT HD ZOOM	ARV120 (GREEN)
CF-H260AZL/I	ADULT HD ZOOM	ARV120 (GREEN)
CF-H260DL/I	ADULT HD SCOPE GUIDE	ARV120 (GREEN)
CF-HQ190L/I	ADULT HD	ARV120 (GREEN)

SCOPE MODEL	OLYMPUS	ENDOCUFF VISION MODEL
CF-HQ290L/I	ADULT HD DUAL FOCUS	ARV120 (GREEN)
CF-Q140L/I	ADULT	ARV120 (GREEN)
CF-Q160AL/I	ADULT	ARV120 (GREEN)
CF-Q160DL/I	ADULT	ARV120 (GREEN)
CF-Q160L/I	ADULT	ARV120 (GREEN)
CF-Q160S	ADULT	ARV120 (GREEN)
CF-Q165L/I	ADULT	ARV120 (GREEN)
CF-Q180 AL/I	ADULT	ARV120 (GREEN)
CF-Q240AL/I	ADULT	ARV120 (GREEN)
CF-Q240L/I	ADULT	ARV120 (GREEN)
CF-Q240ZL/I	ADULT	ARV120 (GREEN)
PCF-160AL/I	PAEDIATRIC	ARV130 (PURPLE)
PCF-H180AL/I	PAEDIATRIC HD	ARV130 (PURPLE)
PCF-H190L/I	PAEDIATRIC	ARV130 (PURPLE)
PCF-H290L/I	PAEDIATRIC HD	ARV130 (PURPLE)
PCF-Q180AL/I	PAEDIATRIC HD	ARV130 (PURPLE)
PCF-Q260AL/I	PAEDIATRIC	ARV130 (PURPLE)
PCF-Q260AZI	PAEDIATRIC ZOOM	ARV130 (PURPLE)
CF-H180DL/I	ADULT HD SCOPE GUIDE	ARV140 (ORANGE)
SCOPE MODEL	PENTAX	ENDOCUFF VISION MODEL
EC38-i10M/M2/F/F2/L/NF	ADULT	ARV110 (BLUE)
EC-3890FI/LI/Fi2(HD)/Mi2	ADULT	ARV110 (BLUE)
EC-3890MZI/FZI/LZI/SI2/TFK	ADULT	ARV110 (BLUE)
EC-3890TLK	DUAL CHANNEL ADULT	ARV110 (BLUE)
EC-3890TMK/MK2/FK/FK2/LK	ADULT	ARV110 (BLUE)
EC-380LKP/FKP/FK2P/MK2P	ADULT	ARV120 (GREEN)
EC-3830LK/P/FK2	ADULT	ARV120 (GREEN)
EC-3840FK2	ADULT	ARV120 (GREEN)
EC3870LK/LZK/TLK/FK2/MK2	ADULT	ARV120 (GREEN)
EC-3880FK2	ADULT	ARV120 (GREEN)
EC-38LV	ADULT	ARV120 (GREEN)
EC-3970LK	ADULT	ARV120 (GREEN)
EC3470LK	PAEDIATRIC	ARV130 (PURPLE)
EC-3490 FI	PAEDIATRIC HD	ARV130 (PURPLE)
EC-3490 MK/FK/LK/LI/TLI/TFI	PAEDIATRIC	ARV130 (PURPLE)
EC-3490FK	PAEDIATRIC HD	ARV130 (PURPLE)
EC-34-i10M/F/L	PAEDIATRIC HD	ARV130 (PURPLE)

NOTE: FOR CANADA ONLY

Not all medical devices listed are licenced in accordance with Canadian law. Please contact your Olympus representative to find out what is registered for sale in your market.

For more information please contact:
Arc Medical Design Ltd, 43 Park Place, Leeds, LS1 2RY, United Kingdom
info@arcmedicaldesign.com

FOR LATEST COMPATIBILITY SEE WWW.ARCMEDICALDESIGN.COM



ENDOCUFF VISION®

Compatibility Schedule Feb 2019

Issue Ref: CSFEB2019

SCOPE MODEL	STORZ	ENDOCUFF VISION MODEL
13924-PKS	ADULT	ARV120 (GREEN)
13925-PKS	ADULT	ARV120 (GREEN)
13906-PKS	PAEDIATRIC	ARV130 (PURPLE)
13907-PKS	PAEDIATRIC	ARV130 (PURPLE)
13904-PKS	ADULT	ARV140 (ORANGE)
13905-PKS	ADULT	ARV140 (ORANGE)

NOTE: FOR CANADA ONLY
Not all medical devices listed are licenced in accordance with Canadian law. Please contact your Olympus representative to find out what is registered for sale in your market.

For more information please contact:
Arc Medical Design Ltd, 43 Park Place, Leeds, LS1 2RY, United Kingdom
info@arcmedicaldesign.com

FOR LATEST COMPATIBILITY SEE WWW.ARCMEDICALDESIGN.COM